



Product Specification **TS/5-10**, 5ml Amies with charcoal transport medium in tube (blue cap) with orange cap cotton minitipped plasticised paper shaft swab in a peel pouch

Product Code	TS/5-10
Description	5ml Amies with charcoal transport medium in tube (blue cap) with orange cap cotton minitipped plasticised paper shaft swab in a peel pouch
Intended Use	Bacteriological samples only. Collection and transport of clinical specimens
Components	Blue cap; Polypropylene Tube; Orange cap; cotton tip plasticised paper shaft swab; Non-toxic glue; Label & film; Amies Transport Medium;
Contents	De-ionised water; Amies Transport Medium;
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	2 years from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.ii) Organism Recovery

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 08.04.2008  
ProductSpec\_TS5-10

