



## Product Specification Sheet

### **TS/5-18A**

4ml Amies with charcoal transport medium in tubewith black cap pure cotton tipped polystyrene shaftswab in a peel pouch

Product Code	TS/5-18A
Description	4ml Amies with charcoal transport medium in tubewith black cap pure cotton tipped polystyrene shaftswab in a peel pouch
Intended Use	Bacteriological samples only. Collection and transport of clinical specimens
Components	Black cap Polypropylene Tube (80 mm) Black cap cotton tipped plastic swab (shaft length80 mm) Non-toxic glue Label & film
Contents	De-ionised waterAmies Charcoal Transport MediumpH – 7.15 ± 0.2
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect fromlight, heat and moisture in tightly closedoriginal container
Pack Size	100 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.ii)Organism Recovery (N.gonorrhoeae,H.influenzae)

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 30.01.2008

ProductSpec\_TS5-18A