



Product Specification Sheet

TS/5-2

5ml Amies with charcoal transport medium in tube(black cap) with black cap cotton wool tipped(wood shaft) transport swab in a peel pouch

Product Code	TS/5-2
Description	5ml Amies with charcoal transport medium in tube(black cap) with black cap cotton wool tipped(wood shaft) transport swab in a peel pouch
Intended Use	This product is intended for bacteriological samples only DO NOT USE THIS PRODUCT WHEN CHLAMYDIAL INFECTION IS EXPECTED
Components	Black cap Polypropylene Tube Black cap cotton/wood swab Non-toxic glue Label & film Transport medium
Contents	De-ionised water Sodium chloride Disodium hydrogen phosphate Sodium thioglycollate Calcium chloride Magnesium chloride Potassium phosphate Agar
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism Recovery

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 31.10.2007

ProductSpec_TS5-2