



Product Specification Sheet

TS/5-20

5ml Amies with charcoal transport medium intube with orange cap cotton tipped twistedwire shaft swab in a peel pouchLatex free

| | |
|-------------------------|--|
| Product Code | TS/5-20 |
| Description | 5ml Amies with charcoal transport medium intube with orange cap cotton tipped twistedwire shaft swab in a peel pouchLatex free |
| Intended Use | Bacteriological samples only. Collection andtransport of clinical specimens. |
| Components | Polypropylene TubeBlack CapOrange cap cotton/twisted wire swabNon-toxic glueLabel & film Transport medium |
| Contents | De-ionised waterAmies Charcoal Transport MediumpH – 7.15 ± 0.2 |
| Method of Sterilisation | Sterilised by gamma irradiation at dosage of25-42KGy |
| Shelf Life | 24 months from date of manufacture |
| Storage Conditions | 4-25 °C; in dry and ventilated area; protect fromlight, heat and moisture in tightly closed original container |
| Pack Size | 100 in a box |
| Other | Quality control testing details on final product:i) Validation of irradiation sterility.ii) Organism Recovery |

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 09.04.2008

ProductSpec_TS5-20