



Product Specification **TS/5-28**, 5ml Amies transport medium with charcoal in tube with blue cap pure viscose tipped (polystyreneshaft) dual swab in a peel pouch

Product Code	TS/5-28
Description	5ml Amies transport medium with charcoal in tube with blue cap pure viscose tipped (polystyreneshaft) dual swab in a peel pouch
Intended Use	Bacteriological samples only. Collection and transport of clinical specimens
Components	Blue cap; Polypropylene Tube; Dual swab: Blue cap viscose tip, polystyrene (breakable) shaft; Non-toxic glue; Label & film;
Contents	De-ionised water; Amies Transport Medium; pH – 7.15 + 0.2
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	2 years from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.ii) Organism Recovery

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 15.02.2008
ProductSpec_TS5-28

