



Product Specification **TS/5-59**, 5ml Marine Salt transport medium in tube with blue standard cap viscose tipped Polystyrene shaft swab in a peel pouch

Product Code	TS/5-59
Description	5ml Marine Salt transport medium in tube with blue standard cap viscose tipped Polystyrene shaft swab in a peel pouch
Intended Use	Bacteriological samples only. Collection and transport of Marine or halophilic organisms.
Components	Polypropylene Tube; Blue Cap x2; Viscose/Polystyrene swab; Non-toxic glue; Label & film; Marine Transport medium
Contents	De-ionised water; Charcoal; Marine Salt Mix;
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	2 years from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.ii) Organism Recovery

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 09.04.2007
ProductSpec_TS5-59

