



Product Specification Sheet

TS/5-6

5ml Amies with charcoal transport medium in tube (black cap). Includes aluminium shaft, cotton mini-tip swab (orange cap) in a peel pouch.

Product Code	TS/5-6
Description	5ml Amies with charcoal transport medium in tube (black cap). Includes aluminium shaft, cotton mini-tip swab (orange cap) in a peel pouch.
Intended Use	Bacteriological samples only. Collection and transport of clinical specimens
Components	Black cap Polypropylene Tube Orange cap cotton/aluminium swab Non-toxic glue Label & film
Contents	De-ionised water Amies Charcoal Transport Medium
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42 KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism Recovery

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 3 Date : 08.02.2009

ProductSpec_TS5-6