



Product Specification Sheet

TS/5-8

5ml Amies with charcoal transport medium in tube (black cap) with orange cap cotton minitipped plasticised paper shaft swab in a peel pouch

Product Code	TS/5-8
Description	5ml Amies with charcoal transport medium in tube (black cap) with orange cap cotton minitipped plasticised paper shaft swab in a peel pouch
Intended Use	Bacteriological samples only. Collection and transport of clinical specimens
Components	Black cap Polypropylene Tube Orange cap cotton tip plasticised paper shaft swab Non-toxic glue Label & film
Contents	De-ionised water Amies Charcoal Transport Medium
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42 KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism Recovery

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 08.04.2008

ProductSpec_TS5-8