



Product Specification **TS/5-47**, Liquid Stuarts Transport Medium in tube with viscose/polystyrene swab in a peel pouch

Product Code	TS/5-47
Description	Liquid Stuarts Transport Medium in tube with viscose/polystyrene swab in a peel pouch
Intended Use	Collection and transport of specimens including N. gonorrhoeae
Components	Blue standard cap; white viscose/polystyrene swab; Non-toxic glue; Tube Label & film; Stuarts Transport medium;
Contents	De-ionised water; Liquid Stuarts Transport Medium 16g; pH – 7.4 ± 0.2
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	2 years from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.ii) Organism Recovery

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 29.03.2007
ProductSpec_TS5-47

