



## Product Specification Sheet

### **M01B004**

Malonate/Phenylalanine Broth, 4ml in Glass Bijou, Terminally sterilised

Product Code	M01B004
Description	Malonate/Phenylalanine Broth, 4ml in Glass Bijou, Terminally sterilised
Intended Use	Use as part of classic testing regime for identification of unknown Microorganisms.
Components	Bijou bottle; Plastic Cap; Polystyrene tray; Label & film;
Contents	Yeast extract; Ammonium sulphate; Dipotassium hydrogen phosphate; Potassium dihydrogen phosphate; Sodium chloride; Sodium malonate; L-Phenylalanine; Bromothymol blue
Method of Sterilisation	Autoclaving at 121°C, 15 min
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 18.09.2012

ProductSpec\_M01B004