



Product Specification Sheet

M02B008

Amies Charcoal Transport Medium, terminally sterilised, 8ml in bijou

Product Code	M02B008
Description	Amies Charcoal Transport Medium, terminally sterilised, 8ml in bijou
Intended Use	Transporting & preserving microbiological samples upto 48hrs
Components	Bijou bottle; Plastic Cap; Polystyrene tray; Label & film;
Contents	Sodium chloride; Potassium chloride; Calcium chloride; Magnesium Chloride; Monopotassium Phosphate; Disodium Phosphate; Sodium Thioglycollate; Charcoal
Method of Sterilisation	Autoclaving at 121°C, 15 min
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product:i)

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 18.09.2012

ProductSpec_M02B008