



Product Specification Sheet

M35B008

Trichomonas Broth, 8ml in bijou

Product Code	M35B008
Description	Trichomonas Broth, 8ml in bijou
Intended Use	Isolating Trichomonas vaginalis. Will also isolate yeasts
Components	Bijou bottle; Plastic Cap; Polystyrene tray; Label & film;
Contents	Tryptone; Peptone; Yeast Extract; Maltose; L-Cysteine; Ascorbic acid; Potassium Phosphate monobasic; Potassium Phosphate bibasic; Chloramphenicol
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 20-46 kGy
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product:i)

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 3 Date : 26.07.2013

ProductSpec_M35B008