



Product Specification **M62B00S**, Nutrient Agar Slopes in Bijou bottles

Product Code	M62B00S
Description	Nutrient Agar Slopes in Bijou bottles
Intended Use	General purpose medium for non-fastidious organisms
Components	Small Bijou bottle, tubular glass Plastic Cap Polystyrene tray Label & film
Contents	Deionised water; Nutrient Agar Powder 28g;
Method of Sterilisation	Autoclaving at 121 °C, 15 min
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25 °C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product: i) Growth of E.coli; Growth of S. aureus; No growth in uninoculated control after 24h at 37 °C iv) No growth in uninoculated control after 7d at 22 °C v) pH = 7.3 + 0.2

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 07.09.2010
ProductSpec_M09B00S

