



Product Specification **M65B004**, Peptone Water Sugars – Glucose, contains Andrade indicator; 4ml in Bijou with Durham tube

Product Code	M65B004
Description	Peptone Water Sugars – Glucose, contains Andrade indicator; 4ml in Bijou with Durham tube
Intended Use	Medium for fermentation studies. Nutrient base with Glucose; contains Andrade's indicator
Components	Bijou bottle; Plastic Cap; Durham tube; Polystyrene tray Label & film
Contents	Peptone Water Sugars - Glucose: Yeast Extract; Tryptone; Meat peptone; Glucose; Andrade indicator; Deionised water; pH – 7.6 + 0.2
Method of Sterilisation	Autoclaving at 121 °C, 15 min
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25 °C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product:i) Positive for E.coli (pink)ii) Negative for P.aeruginosa (no colourchange)iii) No growth in uninoculated control after 24h at 37 °Civ) No growth in uninoculated control after 7d at 22 °C Appearance - Clear pale pink liquid

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 07.09.2010
ProductSpec_M65B064

