



Product Specification **R50-05**, Potassium Tellurite (3.5%), sterile 2ml vials. Pack of 10

Product Code	R50-05
Description	Potassium Tellurite (3.5%), sterile 2ml vials. Pack of 10
Intended Use	Reagent for Bacterial selection and identification, especially useful for <i>S. aureus</i>
Components	Metal crimped cap; Rubber stopper; Glass vial; Potassium Tellurite
Contents	Potassium Tellurite Solution
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42K Gy
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	10 vials in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.ii) Organism Recovery

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 15.02.2007

ProductSpec_R50-05

