



Product Specification **R50-07**, Kovacs Reagent, 4 x 25mls in dropper bottles.

Product Code	R50-07
Description	Kovacs Reagent, 4 x 25mls in dropper bottles.
Intended Use	Reagent for identification, especially useful for E.coli
Components	Plastic cap with rubber teat; Brown Glass Bottle; Kovacs Reagent; Minigrip bag
Contents	HCL; DMAB; Butanol;
Method of Sterilisation	Not applicable
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	4 x 25mls
Other	Quality control testing details on final product:i) Validation of irradiation sterility.ii) Organism Recovery

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 15.02.2007

ProductSpec\_R50-07

