



Product Specification **TS/15-F**, Sterile 500 ml Empty Sampling Container

|                         |                                                                                                                       |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Code            | TS/15-F                                                                                                               |
| Description             | Sterile 500 ml Empty Sampling Container                                                                               |
| Intended Use            | Sampling of potable and swimming pool water samples for microbiological analysis or for general purpose applications. |
| Components              | HDPE Bottle HDPE Cap Label                                                                                            |
| Contents                | Not applicable as supplied empty.                                                                                     |
| Method of Sterilisation | Sterilized by gamma irradiation at dosage of 25-42 kGy                                                                |
| Shelf Life              | 2 years from the date of manufacture                                                                                  |
| Storage Conditions      | 4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container       |
| Pack Size               | 100 in a box                                                                                                          |
| Other                   | Quality control testing details on final product: 1)<br>Validation of irradiation sterility                           |

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 07.05.2009

ProductSpec\_TS15-F

