



Product Specification Sheet

TS/15-Z

Sterile 1000 ml Water Sampling Container with Sodium Thiosulfate (100 mg)

Product Code	TS/15-Z
Description	Sterile 1000 ml Water Sampling Container with Sodium Thiosulfate (100 mg)
Intended Use	Sampling of potable and swimming pool watersamples for microbiological analysis. Sodium thiosulphate neutralises chlorine (up to 28mg/L)
Components	HDPE Bottle HDPE Cap Label
Contents	Sodium thiosulphate- 100 mg (1ml of 10%Na ₂ S ₂ O ₃)
Method of Sterilisation	Sterilized by gamma irradiation at dosage of 25-42 kGy
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25°C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product:1) Validation of irradiation sterility

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 07.05.2009

ProductSpec_TS15-Z