



Product Specification Sheet

TS/163-B

Sterile 1000 ml Empty Sampling Container

| | |
|-------------------------|--|
| Product Code | TS/163-B |
| Description | Sterile 1000 ml Empty Sampling Container |
| Intended Use | Multipurpose Use |
| Components | HDPE Bottle HDPE Cap Label |
| Contents | Not applicable as supplied empty. |
| Method of Sterilisation | Sterilized by gamma irradiation at dosage of 25-42 kGy |
| Shelf Life | 24 months from date of manufacture |
| Storage Conditions | 4-25°C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container |
| Pack Size | 100 in a box |
| Other | Quality control testing details on final product: 1) Validation of irradiation sterility |

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 14.09.2012

ProductSpec_TS163-B