



## Product Specification Sheet

### **TS/19-3DNS**

3" Polyester tip /Polystyrene shaft swab.Single in the tube. Non-sterile.

Product Code	TS/19-3DNS
Description	3" Polyester tip /Polystyrene shaft swab.Single in the tube. Non-sterile.
Intended Use	General swabbing or disgnostic applications where sterility is not essential.
Components	Blue standard cap 3" Polyester tip/polystyrene shaft swab Polypropylene Tube Box Label & film
Contents	Not applicable.
Method of Sterilisation	Non-sterile
Shelf Life	N/A
Storage Conditions	Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product:In-process inspection only

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 08.08.2007

ProductSpec\_TS19-3DNS