



## Product Specification Sheet

### **TS/164-F**

Sterile 2500 ml Empty Sampling Container

Product Code	TS/164-F
Description	Sterile 2500 ml Empty Sampling Container
Intended Use	Sampling of potable and swimming pool water samples for microbiological analysis or for general purpose applications.
Components	HDPE Bottle HDPE Cap Label
Contents	Not applicable as supplied empty.
Method of Sterilisation	Sterilized by gamma irradiation at dosage of 25-42 kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	50 in a box
Other	Quality control testing details on final product: 1) Validation of irradiation sterility

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 07.05.2009

ProductSpec\_TS164-F