



Product Specification Sheet

TS/11

Sterile Orange sticks (wooden swab shaft alone) no tip, 100 in sealed Minigrip bag

Product Code	TS/11
Description	Sterile Orange sticks (wooden swab shaft alone) no tip, 100 in sealed Minigrip bag
Intended Use	Sample throughout. General laboratory use only
Components	Wooden shaft swab; Minigrip bag; Label & film;
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42K Gy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	5000 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 09.04.2008

ProductSpec_TS11