



## Product Specification Sheet

### **TS/6-S**

Sterile 3" sharpened tip cotton swab with a wooden shaft in 6" polypropylene outer tube

Product Code	TS/6-S
Description	Sterile 3" sharpened tip cotton swab with a wooden shaft in 6" polypropylene outer tube
Intended Use	Culture swab; general laboratory use; forensic or precision work
Components	6" Polypropylene Tube 3" Swab with Blue cap Cotton tip/wooden shaft Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	500 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 30.05.2008

ProductSpec\_TS6-S