



## Product Specification Sheet

### **TS/2-A**

10 µl Green Culture Loops (firm), 20 in thesealed minigrip bag

Product Code	TS/2-A
Description	10 µl Green Culture Loops (firm), 20 in thesealed minigrip bag
Intended Use	Sample throughout. General laboratory use only
Components	10 µl Plastic Green Loops (firm) Minigrip bag
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect fromlight, heat and moisture in tightly closedoriginal container
Pack Size	1000 in a box
Other	Quality control testing details on final product:Validation of irradiation sterility

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 19.12.2012

ProductSpec\_TS2-A