



## Product Specification Sheet

### **TS/4-D**

1 µl Yellow Culture Loops (soft), 20 in thesealed minigrip bag

Product Code	TS/4-D
Description	1 µl Yellow Culture Loops (soft), 20 in thesealed minigrip bag
Intended Use	Sample throughout. General laboratory use only
Components	1 µl Plastic Yellow Loops (soft) Minigrip bag
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect fromlight, heat and moisture in tightly closedoriginal container
Pack Size	1000 in a box
Other	Quality control testing details on final product:Validation of irradiation sterility

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 19.12.2012

ProductSpec\_TS4-D