



Product Specification Sheet

TS/5-B

10 µl Blue Culture Loops (soft), 50 in a sealed minigrip bag

Product Code	TS/5-B
Description	10 µl Blue Culture Loops (soft), 50 in a sealed minigrip bag
Intended Use	Sample throughout. General laboratory use only
Components	10 µl Plastic Blue Loops (soft) Minigrip bag
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 12.5-42K Gy
Shelf Life	36 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	1000 in a box
Other	Quality control testing details on final product: Validation of irradiation sterility

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 19.12.2012

ProductSpec_TS5-B