



## Product Specification Sheet

### **TS/30-ES**

Blue L shaped Spreaders (soft), 25 in a sealed minigrip bag

Product Code	TS/30-ES
Description	Blue L shaped Spreaders (soft), 25 in a sealed minigrip bag
Intended Use	Sample throughout or Protect. General laboratory use only
Components	L shaped Blue Plastic (soft) Minigrip bag
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42K Gy
Shelf Life	36 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	1000 in a box
Other	Quality control testing details on final product: Validation of irradiation sterility

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 19.12.2012

ProductSpec\_TS30-ES